

Insulin Syringes

- Boxes of 100, 10 sterile poly bags of 10 syringes for your convenience.
- Clear bold gradation for easy reading.
- Needles are UV bonded to ensure stability.
- Permanent barrel markings assure accurate dosing injection.
- Lubricated needle for smoother injection.
- Easy to read tri-lingual packaging.

U-40

| | Item Number | Product Description | UPC | Box Quantity | Case Quantity |
|--------------------|-------------|---------------------|--------------|--------------|---------------|
| Half Unit Markings | 12-4903 | 29G x 1/2" 3/10cc | 786227781051 | 100 ct. | 5/bxs |
| | 12-4905 | 29G x 1/2" 1/2cc | 786227780559 | 100 ct. | 5/bxs |
| | 12-4910 | 29G x 1/2" 1cc | 786227790558 | 100 ct. | 5/bxs |

U-100

| | Item Number | Product Description | UPC | Box Quantity | Case Quantity |
|--------------------|-------------|---------------------|--------------|--------------|---------------|
| | 12-7805 | 28G x 1/2" 1/2cc | 786227780559 | 100 ct. | 5/bxs |
| | 12-7810 | 28G x 1/2" 1cc | 786227781051 | 100 ct. | 5/bxs |
| | 12-7903 | 29G x 1/2" 3/10cc | 786227790350 | 100 ct. | 5/bxs |
| | 12-7905 | 29G x 1/2" 1/2cc | 786227790558 | 100 ct. | 5/bxs |
| | 12-7910 | 29G x 1/2" 1cc | 786227791050 | 100 ct. | 5/bxs |
| Half Unit Markings | 12-7653 | 31G x 5/16" 3/10cc | 786227765358 | 100 ct. | 5/bxs |
| | 12-7655 | 31G x 5/16" 1/2cc | 786227765556 | 100 ct. | 5/bxs |
| | 12-7657 | 31G x 5/16" 1cc | 786227765754 | 100 ct. | 5/bxs |

Pen Needles

- Lubricated needle provides the best possible injection comfort.
- Compatible with all major brand insulin delivery devices.
- Easy to read tri-lingual packaging.

| | Item Number | Product Description | UPC | Box Quantity | Case Quantity |
|--|-------------|---------------------|---------------|--------------|---------------|
| | 12-5905 | 29G x 1/2" | 7862277590554 | 100 ct. | 12/bxs |

