

Safety Lancets

- Easy to use, pressure activated
- Designed for both patients and healthcare professionals
- Needle retracts after use to leave no possibility of reuse
- New easy to read, bi-lingual packaging

Dual use 18G blade lancet is ideal for tests requiring larger blood samples. May be used on finger or heel.

High flow 21G lancet is used when a larger blood sample is needed. Used for testing blood gas, cholesterol and hemoglobin.

Medium flow 23G lancet is perfect for use when testing blood glucose.

Low flow 30G offers minimal pain as the thinnest lancet in the line.

18 Gauge *Dual Use* High Flow- Finger stick
Low Flow- Heel stick

21 Gauge High Flow

23 Gauge Medium Flow

30 Gauge Low Flow

1.8 mm Depth Setting

2.2 mm Depth Setting

2.2 mm Depth Setting

1.6 mm Depth Setting



18G x 1.8mm
Item # 28-1810
NDC # 86227-0181-05

21G x 2.2mm
Item # 28-2110
NDC # 86227-0211-05

23G x 2.2mm
Item # 28-2310
NDC # 86227-0231-05

30G x 1.6mm
Item # 28-3010
NDC # 86227-0301-15

