

# SureComfort™ Our Diabetic Product Line

## is defined by **Comfort, Quality** and **Affordability**



### Insulin Syringes

SureComfort lubricated needles promote injection comfort.

Our quality syringes have:

- Clear, bold gradations for reading clarity helps patients and care givers deliver an accurate dose.
- Oversized finger grip supports better injection control.
- Each shelf box contains one hundred sterile, non-toxic, pyrogen-free insulin syringes packaged in poly bags of ten.
- Not made with natural rubber Latex.
- Integrated needles minimize drug loss.
- Needles are UV bonded to ensure stability.
- Listed in First Data Bank.

	Item Number	Insulin Syringe Description	NDC	Box Quantity	Case Quantity	Box Dimensions	Weight lb.	Case Dimensions	Weight lb.	Master Case Dimensions	Weight lb.	Pallet Dimensions
28G	22-8005	28G x 1/2" 1/2cc	86227 0800 55	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-8010	28G x 1/2" 1cc	86227 0801 05	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
29G	22-9003	29G x 1/2" 3/10cc	86227 0900 35	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-9005	29G x 1/2" 1/2cc	86227 0900 55	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-9010	29G x 1/2" 1cc	86227 0901 05	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
30G	22-6003	30G x 5/16" 3/10cc	86227 0700 35	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-6005	30G x 5/16" 1/2cc	86227 0700 55	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
30G	22-6010	30G x 5/16" 1cc	86227 0701 05	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-6203	30G x 1/2" 3/10cc	86227 0620 35	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-6205	30G x 1/2" 1/2cc	86227 0620 55	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
31G	22-6210	30G x 1/2" 1cc	86227 0621 05	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-6403	31G x 1/4" 3/10cc	86227 6403 05	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-6405	31G x 1/4" 1/2cc	86227 6405 05	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
31G	22-6410	31G x 1/4" 1cc	86227 6410 05	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-6503	31G x 5/16" 3/10cc	86227 0650 35	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-6504	31G x 5/16" 3/10cc*	86227 0650 45	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-6505	31G x 5/16" 1/2cc	86227 0650 55	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85
	22-6510	31G x 5/16" 1cc	86227 0651 05	100 ct.	5/bxs	5.5 x 7 x 3	0.95	7 x 17.15 x 6	5.25	21.65 x 17.15 x 12.75	31.6	48 x 40 x 85

\*HUM - Half Unit Markings



### Pen Needles

SureComfort pen needles provide optimum comfort and convenience. Designed to be compatible with all insulin pens available in the U.S.

- Five comfortable choices for insulin pen users.
- Boxes contain one hundred sterile, non-toxic, pyrogen free pen needles.
- Not made with natural rubber Latex.

	Item Number	Pen Needle Description	NDC	Box Quantity	Case Quantity	Box Dimensions	Weight lb.	Case Dimensions	Weight lb.	Master Case Dimensions	Weight lb.	Pallet Dimensions
29G	24-1010	29G x 1/2" 12mm	86227 0101 05	100 ct.	12/bxs	5 x 5 x 3	0.25	5 x 15.75 x 6	3.55	18 x 15.75 x 12.75	24	48 x 40 x 85
30G	24-1115	30G x 5/16" 8mm	86227 0111 55	100 ct.	12/bxs	5 x 5 x 3	0.25	5 x 15.75 x 6	3.55	18 x 15.75 x 12.75	24	48 x 40 x 85
31G	24-1215	31G x 5/16" 8mm	86227 0121 55	100 ct.	12/bxs	5 x 5 x 3	0.25	5 x 15.75 x 6	3.55	18 x 15.75 x 12.75	24	48 x 40 x 85
	24-1210	31G x 3/16" 5mm	86227 0121 05	100 ct.	12/bxs	5 x 5 x 3	0.25	5 x 15.75 x 6	3.55	18 x 15.75 x 12.75	24	48 x 40 x 85
32G	24-1305	32G x 1/4" 6mm	86227 0130 55	100 ct.	12/bxs	5 x 5 x 3	0.25	5 x 15.75 x 6	3.55	18 x 15.75 x 12.75	24	48 x 40 x 85
	24-1308	32G x 5/32" 4mm	86227 0130 85	100 ct.	12/bxs	5 x 5 x 3	0.25	5 x 15.75 x 6	3.55	18 x 15.75 x 12.75	24	48 x 40 x 85



### Safety Pen Needles

This Passive/Automatic safety shield system is designed with a red indicator that designates the safety pen needle is ready to be used. The audible and tactile "click" activated upon injection indicates the safety feature is ready to activate without further operation avoiding additional contact.

	Item Number	Safety Pen Needle Description	NDC	Box Quantity	Case Quantity	Box Dimensions	Weight lb.	Case Dimensions	Weight lb.	Master Case Dimensions	Weight lb.	Pallet Dimensions
31G	74-3106	31G x 1/4" 6mm	86227 7402 05	100 ct.	10/bxs	5.6 x 3.5 x 5	1	18 x 6 x 11	11	19 x 13 x 19	48	48 x 40 x 85
32G	74-3204	32G x 3/32" 4mm	86227 7401 05	100 ct.	10/bxs	5.6 x 3.5 x 5	1	18 x 6 x 11	11	19 x 13 x 19	48	48 x 40 x 85



### Lancets

Find maximum comfort in a tri-beveled SureComfort lancet. Designed for universal lancing devices, these lancets are easy to handle and independently sealed to maintain sterility.

	Item Number	Lancet Description	NDC	Box Quantity	Case Quantity	Box Dimensions	Weight lb.	Case Dimensions	Weight lb.	Master Case Dimensions	Weight lb.	Pallet Dimensions
28G	26-2810	28G .36mm	86227 0281 05	100 ct.	50/bxs	2.4 x 1.6 x 3.1	0.1	-	-	12.8 x 8.8 x 7.6	7.7	48 x 40 x 85
30G	26-3010	30G .30mm	86227 0301 05	100 ct.	50/bxs	2.4 x 1.6 x 3.1	0.1	-	-	12.8 x 8.8 x 7.6	7.7	48 x 40 x 85



### Safety Lancets

SureComfort's pressure activated, easy to use safety lancets are designed for both patients and healthcare professionals. Needle retracts after use minimizing the possibility of reuse and avoiding additional contact.

	Item Number	Safety Lancets Description	NDC	Box Quantity	Case Quantity	Box Dimensions	Weight lb.	Case Dimensions	Weight lb.	Master Case Dimensions	Weight lb.	Pallet Dimensions
18G	28-1810	18G x 1.8mm	86227 0181 05	100 ct.	10/bxs	4.7 x 3.9 x 2.4	0.6	3 x 5 x 10	7	20.3 x 13.5 x 10.2	26.5	48 x 40 x 85
21G	28-2110	21G x 2.2mm	86227 0211 05	100 ct.	10/bxs	4.7 x 3.9 x 2.4	0.6	3 x 5 x 10	7	20.3 x 13.5 x 10.2	26.5	48 x 40 x 85
23G	28-2310	23G x 2.2mm	86227 0231 05	100 ct.	10/bxs	4.7 x 3.9 x 2.4	0.6	3 x 5 x 10	7	20.3 x 13.5 x 10.2	26.5	48 x 40 x 85
30G	28-3010	30G x 1.6mm	86227 0301 15	100 ct.	10/bxs	4.7 x 3.9 x 2.4	0.6	3 x 5 x 10	7	20.3 x 13.5 x 10.2	26.5	48 x 40 x 85



### Lancing Pen Device

SureComfort's push-and-click style lancing pen features adjustable settings with an audible and tactile "click" indicating activation. Box contains one lancing pen with easy-to-follow instructions.

	Item Number	Lancing Pen Description	NDC	Box Quantity	Case Quantity	Box Dimensions	Weight lb.	Case Dimensions	Weight lb.	Master Case Dimensions	Weight lb.	Pallet Dimensions
1ea	26-5220	Push & Click	86227 0522 05	1 ct.	25/bxs	1.1 x 1.1 x 5.25	0.6	5.5 x 5.5 x 5.25	1	21 x 12 x 12.5	23	48 x 40 x 85



### Alcohol Prep Pads

SureComfort's moist, four layer thick sterile pads are saturated with 70% v/v isopropyl alcohol. Packaged in air-tight, sealed wrappers for easy storage and convenience.

	Item Number	Alcohol Prep Pads Description	NDC	Box Quantity	Case Quantity	Box Dimensions	Weight lb.	Case Dimensions	Weight lb.	Master Case Dimensions	Weight lb.	Pallet Dimensions
100	26-1260	Single & Sterile	86227 000 00	100 ct.	50/bxs	4 x 2 x 2	0.2	-	-	11.2 x 8.9 x 11.2	14.3	48 x 40 x 85