



# Trusted Quality

ISO 7886-1  
ISO 10993

Country of Origin:  
Korea

- Sterile needles are sealed in individual packages for assurance & convenience.
- Safety needles are available with both Luer Lock and Luer Slip configurations.
- Easy, three-way safety wing activation: forefinger, thumb or hard surface.
- Lubricated needle provides smooth injection.
- Needles are UV bonded to ensure stability.

## 25G & 23G (Variety of Lengths)

## Safety Needles

Item Number	Syringe Description	NDC	Box Quantity	Master Case Quantity	Box Dimensions	Weight lb.	Master Case Dimensions	Weight lb.	Pallet Dimensions
32-8700	25G x 5/8"	86227081005	50 ct.	16/bxs	3.5 x 7 x 3.5	.5	14.75 x 14.25 x 7.5	7.5	48 x 40 x 85
32-8701	25G x 1"	86227081015	50 ct.	16/bxs	4 x 7 x 4	.5	16.25 x 14.25 x 8.5	7.5	48 x 40 x 85
32-8702	25G x 1 1/2"	86227081025	50 ct.	16/bxs	4 x 7 x 4	.5	16.25 x 14.25 x 8.5	7.5	48 x 40 x 85
32-8703	23G x 1"	86227081035	50 ct.	16/bxs	4 x 7 x 4	.5	16.25 x 14.25 x 8.5	7.5	48 x 40 x 85
32-8704	23G x 1 1/2"	86227081045	50 ct.	16/bxs	4 x 7 x 4	.5	16.25 x 14.25 x 8.5	7.5	48 x 40 x 85

## 1cc & 3cc - LUER LOCK SYRINGES

## Conventional Syringes & Safety Needle Combinations

Item Number	Syringe Description	NDC	Box Quantity	Master Case Quantity	Box Dimensions	Weight lb.	Master Case Dimensions	Weight lb.	Pallet Dimensions
35-8200	25G x 5/8" 3cc	86227082005	50 ct.	8/bxs	6 x 1.5 x 4	1	12 x 15.5 x 13	9	48 x 40 x 85
35-8201	25G x 1" 3cc	86227082015	50 ct.	8/bxs	6 x 1.5 x 4	1	12 x 15.5 x 13	9	48 x 40 x 85
35-8202	25G x 1 1/2" 3cc	86227082025	50 ct.	8/bxs	6 x 1.5 x 4	1	12 x 15.5 x 13	9	48 x 40 x 85
35-8203	23G x 1" 1cc	86227082035	50 ct.	12/bxs	7.5 x 9 x 4.5	.8	15.75 x 18 x 9.5	10	48 x 40 x 85
35-8204	23G x 1 1/2" 1cc	86227082045	50 ct.	12/bxs	7.5 x 9 x 4.5	.8	15.75 x 18 x 9.5	10	48 x 40 x 85

We strive to increase your Trust with Quality & Value behind each product produced.