

Personal Protective EQUIPMENT



Disposable Medical Examination Gloves

Designed with a revolutionary new formula of vinyl-based hybrid material containing nitrile. Nitrile uses supple Latex-like elasticity properties and strong bonding which is perfect for most common utilizations.

Nitrile disposable gloves are an indispensable precautionary barrier for your everyday sanitary needs. With a perfect fit, they will cling to your hands like a second skin. Both powder-free and not manufactured with natural rubber Latex, they are the perfect solution for allergy sufferers. Not only are these gloves designed to work with chemicals, they can also be used in the kitchen as they are food-safe.

Features:

- DOH, DEHP-free
- Not made with natural rubber Latex
- Low allergy potential
- Comfortable conforming Latex-like fit
- Effective barrier and extra strength
- Enhanced tactile sensitivity
- Economical protection
- Resistance to acid, alkali and oil

Quality Standards:

- Complies with ASTM D6319
- Complies with EN 455
- Complies with EN 374
- FDA Class 1 registered
- Tensile Strength: 18MPa
- Palm Thickness: 0.08 mm minimum
- Finger Thickness: 0.09 mm minimum
- Weight: 5.0 gm (medium size)
- Length: 230 mm minimum

Technical Information:

- | | |
|--------------------------------|---|
| Type | • Powder-free, non-sterile and not manufactured with natural rubber Latex |
| Model | • Ambidextrous, rolled rim |
| Material | • Compound of Vinyl-based hybrid properties |
| Storage | • Protect from heat, humidity strong light and ozone |
| Impermeability as per EN 455-1 | • AQL 1.5 |

Specifications:

PPE Item Number	PPE Size	PPE Product Description	PPE Box Quantity/UPC	PPE Case Quantity/UPC	PPE Case Dimensions / Weight
GL-503SM	Small	Powder-free Hybrid Gloves	100 / 855303 008501	1,000 / 855303 008546	12.5 x 10.25 x 9.75 13lbs.
GL-503MD	Medium	Powder-free Hybrid Gloves	100 / 855303 008518	1,000 / 855303 008553	12.5 x 10.25 x 9.75 13lbs.
GL-503LG	Large	Powder-free Hybrid Gloves	100 / 855303 008525	1,000 / 855303 008560	12.5 x 10.25 x 9.75 13lbs.
GL-503XL	Extra Large	Powder-free Hybrid Gloves	100 / 855303 008532	1,000 / 855303 008577	12.5 x 10.25 x 9.75 13lbs.



Allison Medical, Inc and its respective designs are trademarks Allison Medical, Inc. Visit Allison Medical's web site at: www.allisonmedical.com

