

# Personal Protective EQUIPMENT



# PPE

## Disposable Medical Mask

Provides a physical barrier to help prevent the direct penetration of droplets and dust. For daily use in medical, industrial, food, beauty, cleaning, household duties and other industries. Non-sterile, single-use and disposable.

### Features:

- 3-ply construction
- Soft and easy to breathe through
- Latex and fiberglass free
- Leak-resistant non-woven fabric
- High-density filter layer
- Elastic ear loops and malleable nose clip
- 98% PFE (particle filtration efficiency)
- 98% BFE (bacterial filtration efficiency)
- Economical protection
- Resistance to acid, alkali and oil

### Quality Standards:

- PE Plastic Nose Wire
- Universal Size

### Compliance Standards:

- Complies with EN 14683:2005

### Specifications:

PPE Item Number	PPE Product Description	PPE Box Quantity/UPC	PPE Case Quantity/UPC	PPE Case Dimensions/Weight
99-0002	Disposable Mask	100ct.	1,600ct.	X - X - X XXlbs.

# PPE

Allison Medical, Inc and its respective designs are trademarks Allison Medical, Inc. Visit Allison Medical's web site at: [www.allisonmedical.com](http://www.allisonmedical.com)

  
**allison** | medical.  
inspired | innovation

608-02-230\_Rev01\_11-22