

SureComfort™

Delivering Quality and Comfort



Safety Pen Needles

	Item Number	Safety Pen Needle Description	NRC Number	Box Quantity	Case Quantity	Box Dimensions & Weight	Case Dimensions & Weight	Master Case Dimensions & Weight
31G	74-3106	31G x 1/4" 6mm	86227 7402 05	100 ct.	10/bxs	5.6 x 3.5 x 5 - 1 lb.	18 x 6 x 11 - 11 lbs.	19 x 13 x 19 - 48 lbs.
32G	74-3204	32G x 3/32" 4mm	86227 7401 05	100 ct.	10/bxs	5.6 x 3.5 x 5 - 1 lb.	18 x 6 x 11 - 11 lbs.	19 x 13 x 19 - 48 lbs.



Pallet sizes are 48x40x85

High Quality in a **PASSIVE/AUTOMATIC SAFETY** Universal Fit* Pen Needle!



EXCLUSIVE FEATURE:

SureComfort's™ "PASSIVE/AUTOMATIC" Safety Pen Needle is a single use system that does not require additional active motion or changing of the hand position for the safety shield to engage after a health-care provider administers the insulin injection. This exclusive feature is designed to minimize accidental needle stick injuries.



NON-CONTACT SYSTEM:

The red indicator designates that the needle is ready to be used. The user will hear an audible and tactile "click" during injection. This indicates the "PASSIVE/AUTOMATIC" safety system has been initiated. The clear safety shield automatically locks in the forward position following an injection without further operation avoiding additional contact.



UNIVERSAL FIT*:

SureComfort™ Safety Pen Needles is the premier quality "PASSIVE/AUTOMATIC" safety pen needle in the market. It is compatible with most leading brand diabetic pens. SureComfort's™ Safety Pen Needle design advances the users control and comfort. SureComfort's™ Safety Pen Needles are available in MICRO and SHORT needle lengths.

Trust and comfort. Every single customer. Every single time.



SureComfort and its respective designs are trademarks of SureComfort, a division of Allison Medical, Inc. Visit Allison Medical's web site at: www.allisonmedical.com

